


Healthometer®  
Professional

Capacity: 400 lb / 180 kg



### Advantages

- ✓ **400 lb / 180 kg Capacity Floor scale**  
The high capacity of this floor scale provides more versatile use
- ✓ **8" Raised Dial Display**  
Large, easy-to-read dial display is raised above the scale base to allow for more footroom on the platform
- ✓ **Weighs in LB/KG**  
Dial display shows weight in both pounds and kilograms providing versatility



0800 436 324  
www.emech.co.nz  
sales@emech.co.nz

**Healthometer®**  
Professional  
www.homscales.com

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l

Capacity: 400 lb / 180 kg



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
160KL	2	17.375" x 15.125" x 21.875" / 3.33	21.5 lb	10892076002364

**Specifications**

- Capacity: 400 lb / 180 kg
- Graduation: 1 lb / 1 kg
- 8" dial (raised 4")
- Platform Size: 11" (w) x 12-1/2" (d) x 3" (h) /  
279 mm (w) x 318 mm (d) x 76 mm (h)
- 1 Year Limited Warranty

emech

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l  
www.homscales.com