

DISPLAYS WEIGHT IN POUNDS, KILOGRAMS & GRAMS

Health o meter[®]
P r o f e s s i o n a l



Advantages

- Motion Detection**
Motion sensing technology compensates for baby's movement and high-resolution readings provide more accurate weight measurements, ideal for the neonatal environment
- Convenient Features**
Plastic tray is easy-to-clean and heavy-duty built-in handles allow for portability
- Built-In Measuring Tape**
Easy-to-read measuring tape allows user to measure both weight and height simultaneously
- Industry Leading 3 Year Warranty**
Health o meter[®] Professional's industry leading warranty protects you against defects of materials or workmanship
- Value-Added Features**
Operating functions (LB/KG/GM Lockout, LB/GM Conversion, LB/KG Conversion, Zero Out/Tare, Reweigh/Recall, Auto Zero, Auto Off) included at no additional cost

Health o meter[®]
P r o f e s s i o n a l
www.hom scales.com

**HIGH RESOLUTION
PORTABLE SCALE**



**NON-CONDUCTIVE
WHEELS**

2210CART

| Model | Shipping Dims. (L x W x H) / Cube | Weight (LB) | UPC Code |
|----------|-----------------------------------|-------------|--------------|
| 2210KL | 27.875" x 18.625" x 7.625" / 2.29 | 30.5 lb | 892076002329 |
| 2210CART | 29" x 21" x 31.625" / 11.15 | 47 lb | |

2210KL



2210KL Specifications

- Capacity: 45 lb / 20 kg / 20,000 gm
- Resolution: 0 – 11300 g (0 – 25 lb) (0 - 10 kg) +/- 1g (+/- 0.1 oz)
> 11300 g (> 25 lb) (> 10 kg) +/- 5 g (+/- 0.1 oz)
- 1 gram accuracy standard
- 7/8" LCD display
- Functions: LB / KG / GM Lockout, LB / GM Conversion, LB / KG Conversion, Zero Out / Tare, Reweigh / Recall, Auto Zero, Auto Off, Motion Detection
- Removable curved clear acrylic tray
- Tray Size: 24" (w) x 15-3/8" (d) x 3" (h) / 610 mm (w) x 391 mm (d) x 76 mm (h)
- Measuring Tape: 0" - 23" / 0 cm - 60 cm
- Heavy duty handles for ease of moving
- 6 C-cell batteries included
- 3 Year Limited Warranty

2210CART Specifications

- Capacity: 200 lb / 90 kg
- Top shelf has handle for easy grip
- Non-conductive wheels / 2 wheels with locks
- Cart Dimensions: 26-1/2" (w) x 18" (d) x 34-1/8" (h) / 673 mm (w) x 457 mm (d) x 867 mm (h)
- Lockable Cabinet: 23" (w) x 16-3/4" (d) x 17" (h) / 584 mm (w) x 427 mm (d) x 432 mm (h)
- Storage Shelf is easy to clean, durable blue plastic
- Use with pediatric scales: 2210KL, 553KL, 549KL
- 1 Year Limited Warranty

