



**High Capacity for  
Bariatric Use**

**Advantages**

- EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter<sup>®</sup> Professional Connectivity Solutions
- “Live” Handrails**  
Advanced technology allows the patient to hold onto handrails for stability while scale accurately calculates weight
- 180° Swivel Display**  
The display head has the ability to rotate /pivot, for the convenience of the administrator, allowing line-of-sight viewing and information privacy
- Low Profile Platform**  
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience
- Value-Added Features**  
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off) included at no additional cost
- Calculates Body Mass Index (BMI)**  
Scale functions include the ability to calculate the patient’s Body Mass Index (BMI)

emech

m e d i c a l

0800 436 324

sales@emech.co.nz

www.emech.co.nz

**Health o meter**<sup>®</sup>

Professional

www.homscales.com



“LIVE” Handrails  
for patient stability

180° Swivel Display



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2101KL	42.125" x 28.625" x 10.5" / 7.33	78 lb	892076002282



**Specifications**

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 26" (w) x 22" (d) x 2 3/8" (h) / 660 mm (w) x 559 mm (d) x 60 mm (h)
- Height Rod: Optional STROD
- Display: 3/4" LCD, 180° Swivel
- Wheels: 2
- Power Source: 110-240V adapter included, 6 AA batteries (not included)
- Connectivity: USB
- Functions: LB / KG Conversion, LB / KG Lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Reweigh, Auto Zero, Auto Off
- 2 Year Limited Warranty