

# Healthometer®

Professional



## Advantages

- Value-Added Features**  
 Operating functions (LB/KG Conversion, Zero Out/Tare, Auto Off, Auto Zero) offered at no additional cost
- Convenient Power Source**  
 Scale is powered by a 6V rechargeable battery and battery charger included with scale
- 1-1/4" LCD Display**  
 The large display shows the user an easy-to-read weight measurement
- Moveable Arm & Foot Rest**  
 The flip-arm and foot rest allows for easier patient access and patient comfort
- Sturdy Construction**  
 The sturdy construction of this scale provides a stable base that does not shake



0800 436 324  
 www.emech.co.nz  
 sales@emech.co.nz

**Healthometer®**  
 Professional  
 www.homscales.com

See reverse side for additional specifications

# Health o meter® P r o f e s s i o n a l



594KL Display

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
594KL	30.125" x 25.375" x 44.625" / 19.74	81.25 lb	892076002688



### Specifications

- Capacity: 600 lb / 270 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size: 18 1/4" (w) x 15" (d) x 17 1/2" (h) / 464 mm (w) x 381 mm (d) x 445 mm (h)
- Display: 1 1/4" LCD
- Wheels: 2
- Power Source: 6 Volt rechargeable battery, charger and battery included
- Functions: LB / KG Conversion, Zero Out / Tare, Auto Zero, Auto Off
- 2 Year Limited Warranty

