

**Health o meter®**  
P r o f e s s i o n a l



**Advantages**

- EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Value-Added Features**  
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- Low Profile Platform**  
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience
- Built-in Height Rod**  
The easy-to-use height rod provides height measurements from 30" to 84 1/4"
- Top-Selling Digital Scale**  
The Health o meter Professional 500KL is the most widely used scale in physician offices nationwide
- Calculates Body Mass Index (BMI)**  
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)

emech

m e d i c a l

0800 436 324

www.emech.co.nz

sales@emech.co.nz

**Health o meter®**  
P r o f e s s i o n a l  
www.homscales.com

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l



| Model | Shipping Dims. (L x W x H) / Cube  | Weight (LB) | UPC Code     |
|-------|------------------------------------|-------------|--------------|
| 500KL | 30.625" x 17.375" x 11.375" / 3.50 | 24.5 lb     | 892076002008 |



**Specifications**

- Capacity: 500 lb / 220 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG Conversion, LB /KG Lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size (w x d x h): 13 3/4" x 16 1/2" x 2 3/8" / 349 mm x 419 mm x 60 mm
- Height Rod: 30" - 84 1/4" / 76 cm - 214 cm
- (2) Wheels for ease of movement
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty

emech  
+ + +  
m e d i c a

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l  
[www.homscales.com](http://www.homscales.com)