

# 5000 SERIES HOMECARE BED



Pin Slot System



Transport Mode (Folded)



Handset



**'IDEAL ALTERNATIVE  
FOR HOME OR  
SHORT-TERM CARE'**

Electric Hi/Low

Electric Backrest

Electric Kneebend

Electric Forward &

Reverse Trendelenburg

SWL 200KG

100mm Braking Castors

Minimum Height - 295mm

Maximum Height - 705mm

P5500 Overall Width - 930mm

Sleeping Surface Width - 2000 x 900mm

P5500W Overall Width - 1070mm

Sleeping Surface Width - 2000 x 1050mm

Overall Length - 2190mm

The Alrick 5000 Pull-A-Part Transportable Series is targeted at easy transportation and storage. Ideal for short term care, the 5000 Series can be assembled and dismantled with ease through the incorporation of a simple slot system which is secured with a pin, and folds compactly for storage. Consequently the 5000 Series is an ideal homecare bed, or as a short-term care alternative.



## ACCESSORIES AVAILABLE TO SUIT 5000 HOMECARE BED



**360HR** Horizontal Side Rails



**380VR** Vertical Side Rails



**400HR** Deluxe Rails



**ASHP/H** Self Help Pole



**AIVP** I.V. Pole



**COT PRO** Padded Side Rail Protectors



**OBT** OverBed Table

## OPTIONS AVAILABLE TO SUIT 5000 HOMECARE BED



**ABB** Battery Backup



**CPR** Emergency Release



**C.EXT5000** 200mm  
Clip On Base Extension



**C.EXT5000W** 200mm  
Clip On Base Extension

## WIDTHS AVAILABLE TO SUIT 5000 HOMECARE BED



**P5500** Standard Single  
900mm Wide



**P5500W** King Single  
1070mm Wide

**ALRICK**

DISTRIBUTOR:

**emech**  
+ + +  
**medical**

0800 436 324 | sales@emech.co.nz