



SOUTHERN CROSS



Supporting you

1300 362 923

www.healthtec.com.au

Universal Examination Table

The SC electric lift Examination Table makes switching to an electric height adjustable table more affordable than ever before.

Offering a generous lifting capacity of 250kg, the motor is covered by an impressive 5 year warranty and moves from HI-LO within 16 seconds enabling the Practitioner to get the patient positioned quickly.

Each end section can be raised to provide a backrest for left and right patient examination and procedures.



Features

- Electric Height Adjustment with footswitch
- X-Frame Lift action minimises pinch and crush points
- Table can be moved similar to moving a wheelbarrow
- Durable powder coated rust resistant frame
- High quality anti-microbial vinyl & CFC free foam
- 5 year warranty on frame & electrics
- Australian designed & manufactured

Specifications

- 250kg Maximum Lifting Capacity (moving)
- 400kg Maximum Working Load (stationary)
- 540mm Minimum Height
- 860mm Maximum Height
- 710mm Table Width
- 1930mm Table Length
- Manual backrest adjustment +70° (both end sections)
- Footswitch Control



The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6

© August 2013




Supporting you

www.healthtec.com.au

PRICE SHEET

Examination Table

	CODE	Qty
Product Description		
SC Universal Examination Table (710 wide)	56221	
Foam & Vinyl Options		
75mm Comfort Foam	6005	
Non-Stock Vinyl Option	6018	
*please refer to our colour chart for vinyl colour options		
Nose Hole Option		
Nose Hole in Backrest	6011-2-7	
Nose Hole Plug	6010	
Armrest Option		
Procedure Arm	5102	
Gynae Option		
Gynae Leg Supports	5114	
Table Protection		
Paper Towel Roll Holder	6017	
Disposable Fitted Bed Sheet – 10 per pack	ANWFIT	
Extended Warranty		
10 years Frame & Electrics	7800	
Freight		
Freight charges are not included in our prices. Please contact your Healthtec representative for freight costs and options.		
Shipping Dimensions	2050mm x 850mm x 650mm	
Shipping Weight	80kg (approximately)	



FRECALL
1300 362 923
 info@healthtec.com.au

Healthtec Pty Ltd ABN 94 224 518 379 TGA AR-00106892-6
 © August 2013

Disclaimer: Due to our policy of constant product improvement, specifications may change without notice. Check with us for the latest specifications and pricing.