



7000 SERIES ULTRA LOW BED



Removable Headboard



Electric Kneebend



Central Locking Castors



'ULTRA LOW 195MM ELECTRIC CARE BED'

Electric Hi/Low

Electric Backrest

Electric Kneebend*
Optional

Vascular Support

SWL 175KG

125mm Braking Castors

Minimum Height - 195mm

Maximum Height - 710mm

Overall Width - 974mm

Overall Length - 2155mm

Sleeping Surface - 2000mm x 880mm

The Alrick 7000 Ultra Low Series caters for individual requirements, while the wide variety of options available makes it one of the most adaptable on the market. The variable height range not only ensures ease of getting in and out of bed, it also minimizes the risk of injury and provides optimum height for care.

*Two models are 7001 with kneebend and vascular support , 7500 without kneebend and vascular support.



*Shown below with Standard Style
Curved Head and Foot Boards*

ACCESSORIES AVAILABLE TO SUIT 7000 ULTRA LOW BED



360HR Horizontal Side Rails



380VR Vertical Side Rails



400HR Deluxe Rails



ASHP Self Help Pole



AIVP I.V. Pole



CLP400 Clamp On Support Rail 400mm



TBR7000 T Bar Buffer



C.EXT 200mm Clip On Extension



BMR Bed Mover Wheel Set



COT PRO Padded Side Rail Protectors



OBT OverBed Table

OPTIONS AVAILABLE TO SUIT 7000 ULTRA LOW BED



ABB Battery Backup



ACM Attendant Control



PO EXT 200mm Pull Out Extension



CLOK 100mm Central Locking Castors



CPR Emergency Release



VSE Electric Vascular Support



UBL Underbed Lighting/
BLH Backlit Handset

HEAD AND FOOT BOARD OPTIONS FOR 7000 ULTRA LOW BED



SS Square Style



QAS Queen Ann Style



SAS Slender Arch Style



ATCS Aluminium Trimmed Crescent Style



DTWS Duet Tone Wooden Style



SMS Sculpted Moulded Style

** Most head and foot board styles are available in a range of colours on request.*

ALRICK

emech
+ + +
m e d i c a l

DISTRIBUTOR:

0800 436 324 | sales@emech.co.nz