

Health o meter[®]
P r o f e s s i o n a l



Advantages

- ✓ **Sanitary Plastic Weighing Surface**
The easy-to-clean plastic tray is ideal for sterile environments
- ✓ **Value-Added Features**
Operating functions (Zero Out, Tare, Auto Zero, Auto Off) included at no additional cost
- ✓ **LB/KG Lockout**
This feature locks the weight measurement unit assuring the weight conversion value is properly recorded when prescribing dosages
- ✓ **Tare Function**
This convenient function allows the user to automatically subtract the weight of blanket or diaper from infant's weight measurement

emech

m e d i c a l

0800 436 324

www.emech.co.nz

sales@emech.co.nz

Health o meter[®]
P r o f e s s i o n a l
www.homscales.com

See reverse side for additional specifications



2210CART
Rolling Scale Cart for 549KL

- Non-conductive wheels / 2 wheels with locks
- Cart Dimensions: 26-1/2" (w) x 18" (d) x 34-1/8" (h)
- Lockable cabinet: 23" (w) x 16-3/4" (d) x 17" (h)
- Storage Shelf
- 1 Year Limited Warranty

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
549KL	1	21.75" x 15.5" x 6" / 1.17	5 lb	892076002336

Specifications

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- 0.6" LCD display
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out / Tare, Auto Zero, Auto Off
- Easy to clean plastic tray
- Tray Size: 19-3/8" (w) x 12-1/8" (d) x 2" (h)
- 9 Volt battery included
- Optional carrying case (order # 64769)
- Optional Rolling Scale Cart (order # 2210CART)
- 2 Year Limited Warranty



Carrying Case for 549KL
Item # 64769