

AARON® 940™

Bovie®

emech
medical
t: 0800 436 324
e: sales@emech.co.nz
w: www.emech.co.nz

ELECTROSURGERY FOR THE OFFICE AND SURGI-CENTER



HIGH-FREQUENCY DESICCATOR

POWER

40 watt high frequency desiccator

FLEXIBILITY

Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar

ACCURACY

More accurate power with Universal Power Control 100-240 VAC

SIMPLICITY

Less to do and less to buy makes this the most cost effective and user friendly unit on the market

COST EFFECTIVE

The Aaron® 940™ is a ground reference unit which means no grounding pad is required

RELIABILITY





The Aaron® 940™ is backed by Bovie's two (2) year manufacturer's limited warranty

AARON® 940™
INCLUDES

3-Button Autoclavable Handpiece	Blunt & Sharp Dermal Tips	Wall Mount Kit
		
A901	A806 & A804 (50 each)	A837

OTHER RECOMMENDED
ACCESSORIES

Footswitch	Reusable Bipolar Forceps*	Forcep Cable
		
A803	A824 & A821	A827EU

Disposable and Reusable Electrodes*			
Ball	Needle	Sterile Sheath	Non-Sterile Sheath
			
ES20	ES02	A910ST	A910

*More styles available

Industry-standard accessories fit most electrosurgical generators

Additional accessories available — contact your Bovie representative for details.

TECHNICAL
SPECIFICATIONS

Technical Specifications...

Dimensions:	Line Voltage:	Line Frequency:
Height 16.0 cm (6.3 in.) Width 22.8 cm (9 in.) Depth 11.4 cm (4.5 in.) Weight < 1.4 kg (< 3 lbs)	Universal Power Supply A940: 100 to 240 VAC	50 - 60 Hz
Output Power:	Output Frequency:	Warranty:
To 40 watts (0-10 watts in 1/10th watt increments on both mono and bipolar)	550 kHz	two (2) year (unit only)

