

AARON<sup>®</sup> 950<sup>™</sup>

*Bovie*<sup>®</sup>

emech  
medical  
t: 0800 436 324  
e: sales@emech.co.nz  
w: www.emech.co.nz

# ELECTROSURGERY FOR THE OFFICE AND SURGERY CENTER



ELECTROSURGICAL GENERATOR

## VERSATILITY

High Frequency Desiccator with 60 watts of Cutting capabilities

## UNIQUE

Two Machines in one convenient package

## SIMPLICITY

Six (6) Doctor presets for fast, easy setup

## ACCURACY

Adjustable in 1/10<sup>th</sup> watt increments in the first 10 watts, for the most accurate output available




## COST EFFECTIVE

No grounding pad required when working in the desiccation mode.

## RELIABILITY

The Aaron 950 is backed by Bovie's two (2) year manufacturer's limited warranty.





AARON® 950™  
INCLUDES

3-Button Autoclavable Handpiece 	Non-Sterile Sheath 	Blunt & Sharp Dermal Tips 
A901	A910 (10 Included)	A806 & A804 (50 each)

Disposable Electrodes 	Solid Disposable Grounding Pad 	Reusable Grounding Cord 	Wall Mount Kit 
ES01 (3 Included), ES02, ES20	A1202 (5 Included)	A1202C	A837

OTHER OPTIONAL  
ACCESSORIES

Footswitch 	Reusable Bipolar Forceps* 	Forcep Cord 
A1203	A825 & A823	A827EU

Disposable and Reusable Electrodes*			Sterile Sheath 
			
Lletz Loop Electrodes*	Resistick II™ Coated Electrodes*	Tungsten Loops*	A910ST

\*More styles and accessories available — contact your Bovie representative for details.  
Industry-standard accessories fit most electrosurgical generators


TECHNICAL  
SPECIFICATIONS

Technical Specifications...

<b>Dimensions:</b>	<b>Line Voltage:</b>	<b>Line Frequency:</b>
Height 16.0 cm (6.3 in.) Width 22.8 cm (9 in.) Depth 11.4 cm (4.5 in.) Weight <3.2 kg (< 7 lbs)	A950: 120 VAC A950-220: 240 VAC	50 - 60 Hz

Monopolar	Output Power	Output Frequency	Rep Rate
Cut	60 W @ 500 Ω	357 kHz ± 50 kHz	N/A
Blend	60 W @ 800 Ω	357 kHz ± 50 kHz	30 kHz ± 5 kHz
Coagulation	60 W @ 1000 Ω	375 kHz ± 50 kHz	60 kHz ± 5 kHz
Fulguration	35 W @ 1000 Ω	575 kHz ± 50 kHz	30 kHz ± 5 kHz
Bipolar	30 W @ 200 Ω	520 kHz ± 14 kHz + 50 kHz	19 kHz ± 5 kHz

**Warranty:** limited two (2) year on unit only



A950 Unit

A808T Tray  
(Optional)

A812 Stand  
(Optional)



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com  
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie® and Aaron® are registered trademarks of Bovie Medical Corporation