


STERILCLAVE® 18/24



Cominox® autoclaves have been projected for facing **binding tasks**: SterilClaves can contain up to **5 big trays** (18 lt.: 180 x 280mm - 24 lt: 180x380mm) or as an alternative up to **3 NORM TYPE cassettes** (310 x 190 x 40mm). All SterilClave models are equipped with “clean” and “dirty” water inspection tank. They are made of resistant **AISI 304 steel**, which ensures the highest level of internal hygiene and durability.

	U.M.	18 / 24 B-BHD	6 / 18 / 24 S Dynamica
Chamber capacity	l	17,15 23,15	4,8 17,15 23,15
Clean water tank capacity	l	2,1 3,1	2,1 2,1 3,1
Used water recycling tank capacity	l	2,1 3,1	-
Water consumption per cycle (minimum - maximum)	cm³	210 - 700 420 - 1190	150 - 250 370 - 590 420 - 1190
Electrical voltage and frequency	V - Hz	230 - 50 (single - phase AC)	230 - 50 (single - phase AC)
Maximum allowed voltage fluctuations	%	±10	±10
Absorbed maximum power	W	2000 2800	1000 2000 2800
Detected noise level in work areas under suitable conditions (medium - peak)	dB(A) - dB(C)	55,7 - 62,7 56 - 63	55,7 - 62,7 55,7 - 62,7 56 - 63
Overall dimensions W x H x D	mm	445 x 390 x 640 445 x 470 x 700	385 x 370 x 535 445 x 390 x 640 445 x 470 x 700
Housing dimensions W x H x D	mm	450 x 395 x 560 450 x 475 x 620	390 x 375 x 505 450 X 395 X 560 450 X 475 X 620
Unit depth with door open	mm	880 940	690 880 940
Usable space dimensions W x H x D	mm	185 x 150 x 285 185 x 150 x 420	185 x 100 x 285 185 x 150 x 285 185 x 150 x 420
Weight (mass)	kg	51 68	38 48,7 62



24 months warranty