



## Discovery

Our best seller over the years the AMW Discovery is the model that suits the most riders. Its large terrain tyres cover rough ground easily and it is a real hill-climber with its 800W motor and 45A batteries. Very popular for community grant applications, this model is priced for complete affordability.

### Features

Motor size:	800 Watt
Battery power:	2 x 45 Amp (12v)
Tyre size:	12"
Dimensions:	1305 x 640mm
Turning radius:	1.6m
Ground clearance:	150mm
Max speed:	10km/hr
Range:	30-35km
Carry capacity:	160kg
Climbing grade:	12° degrees

Brakes:	Electronic regenerative
Tyres:	Pneumatic
Suspension:	Full
Seat:	Adjustable front/back, swivel, height, backrest.
Quick Release:	Dash mounted brake release
Controller:	P&G
Package Includes:	Front/back baskets, rain cover, manual, walking stick holder, walking stick, charger, paint repair, toolkit.

Available in:

