



duplex®

The duplex® stethoscopes from **Riester** combine first-class quality at an attractive price. They offer outstanding acoustic performance for the daily challenges in office practices and hospitals.

The acoustic performance has been considerably improved.

- Double-headed chest-piece available in two types: chromium-plated brass or specially lightweight aluminium, galvanized.
- Special membrane for perfect auscultation, Ø 48 mm (1.89 in), with particularly flat and soft non-chill rims for improved adaptation on the skin.
- Bell, Ø 36 mm (1.42 in)
- Overall length: 77 cm (30.32 in)
- duplex® is supplied with a pair of replacement ear-tips and a replacement membrane.

Order information duplex®:

Aluminium chest-piece

- No. 4001-01 ■ black
No. 4001-02 ■ slate grey
No. 4031-03 ■ blue
No. 4033-05 ■ green

Chest-piece made of chromium-plated brass

- No. 4011-01 ■ black
No. 4011-02 ■ slate grey
No. 4011-03 ■ blue
No. 4011-05 ■ green

duplex®

PRODUCT CHARACTERISTICS