

# P r o f e s s i o n a l

1-800-815-6615 / HOMProCS@homscases.com

**ADULTHR** Height Measuring Rod  
Free-standing stadiometer to measure patient height

**BABYHR** Measuring Rod for Baby  
Wall-mounted stadiometer to measure baby's length

**Range of measurement:** 85 – 210 cm (33 1/2" – 82 5/8")  
**Division of measurement:** 1 mm / 1/16"

**Range of measurement:** 7 – 100 cm / 3 – 39"  
**Division of measurement:** 1 mm / 1/16"

**BABYSLING** Hanging Pediatric Scale  
Portable dial scale with hanging sling

**BMICHARTKG** BMI Measuring Wheel  
Calculates Body Mass Index for children and adults

**Range of measurement:** 15 kg – 25 kg  
**Division of measurement:** 100g

**Product size:** 10 cm diameter  
**Uses kilograms only**