

Choose your **GO₂** Accessories



GO2CC
Black neoprene carrying case with belt loop



GO2L
Black 20" Lanyard



GO2R
Retractable Clip Holder



GO₂
Fingertip Pulse Oximeter

Three colors to Choose from



A Message to Oxygen Users

With Knowledge Comes Peace of Mind

Having the necessary information is essential to successfully managing your oxygen to stay within your clinician-prescribed oxygen saturation and heart rate range. With **GO₂**, you're in control:

- Adjust your oxygen to meet your daily activity levels.
- Get instant feedback on the effectiveness of your breathing techniques.
- Improve communication with your clinician using trusted Nonin Technology.

"I won't leave home without it."
—Edna Fiore, COPD Patient/Advocate

For oxygen patients who love to travel, the new FAA rules allow the use of Portable Oxygen Concentrators (POCs) on commercial airlines.

Go Further! Independence Restored.

A chronic lung disease diagnosis often requires that adjustments be made in your daily activities. However, it does not have to mean a life of limitation! Whether you're looking for peace of mind, to go out for a night at the movies, or to up your activity level by walking those extra blocks around the neighborhood — the **GO₂** fingertip oximeter is the perfect companion! Wherever you decide to go, **GO₂** tucks away in your purse or pocket and provides the assurance that you are within your prescribed oxygen saturation and heart rate range — anywhere, anytime.

"An oximeter gives the patient immediate feedback, a lot of reassurance — reduces anxiety."
—Patty Koff, RRT

Specifications

Oxygen Saturation Display Range	0% to 100
Pulse Rate Display Range	18 to 321 pulses per minute
Oxygen Saturation Declared Accuracy Range (Arms*)	70-100% SpO ₂ ± 2 digits
Low Perfusion Oxygen Saturation Declared Accuracy Range (Arms*)	70-100% SpO ₂ ± 2 digits
Pulse Rate Declared Accuracy Range(Arms*)	20-250 BPM ± 3 digits
Low Perfusion Pulse Rate Declared Accuracy Range (Arms*)	40-240 ± 3 digits
*± 1 Arms represents approximately 68% of measurements.	

Temperature (Operating)	+41° to +104°F (5°C to +40°C)
Storage/Transportation	-22° to +158°F (-30° to +70°C)
Device temperature will not exceed 41°C as measured in a controlled environment.	
Humidity (Operating)	10% to 90% relative humidity, non-condensing
Storage/Transportation	10% to 95% relative humidity, non-condensing
Operating Altitude	Up to 40,000 feet / 12,192 meters
Battery Life (Continuous):	Approximately 2,400 spot checks based on ~21 hours of operation using one AAA-size alkaline battery, calculated at 30 seconds per use.
Battery Life (Storage)	6 months minimum

2 year warranty.
Backlit LCD for low light viewing.
Intended for Home Health Care Only.

Specifications are subject to change without notice.



Nonin Medical, Inc.

13700 1st Avenue North • Plymouth MN • 55441-5443 • USA
Phone: +1.763.553.9968 • Toll Free: 877.577.2635 • Fax: +1.763.577.5521
E-mail: info@nonin.com




www.go2nonin.com
(U.S. Only)



www.nonin.com
(Outside of U.S.)



Authorized Distributor:



©2009 Nonin Medical, Inc.
All trademarks are the property of Nonin Medical, Inc unless otherwise noted.
6960-000-01



Take Control!
With **GO₂** — You Can!



GO₂
Fingertip Pulse Oximeter

Empower Yourself

Nonin Medical, the inventor of fingertip pulse oximetry, introduces the **GO₂** fingertip pulse oximeter. The affordable **GO₂** empowers you with the knowledge of reliable oxygen saturation and heart rate values without compromising the quality and accuracy that have made Nonin the most trusted name in fingertip pulse oximetry.



GO₂

Go Ahead and Get Active!

A comprehensive fitness program is essential to maintaining lung health and quality of life if you have been diagnosed with a chronic lung disease. The highly durable **GO₂** fingertip oximeter is the perfect oxygen saturation and heart rate monitoring device to complement your clinician-devised exercise program. **Know when to increase your liter flow in an instant!** Compact and easy-to-use, simply place **GO₂** on your finger for accurate oxygen saturation and heart rate feedback. Go further with **GO₂**!



“You are not limited as long as you know what your saturation is.”

—Lynn Cole, COPD Patient/Advocate

Quality and Support DO Matter

With so many throwaway products available these days, it’s important to choose a product that not only is consistently accurate, but is also built to last. At Nonin, we are dedicated to providing the highest quality products available at a great value.

Superior Service and Support

Nonin will be here long after the sale and our knowledgeable customer service representatives are always available to answer your questions.

“I relate oximetry to diabetes testing. You wouldn’t want to guess at your sugar levels... You wouldn’t want to guess that maybe three liters of oxygen works when you really need five.”

—Keene Jorgenson, COPD Patient/Advocate

“Patients should all have oximeters to guide their oxygen use during all activities of daily living.”

—Dr. Thomas L. Petty, MD

A Message to Clinicians

Trusted Accuracy

Accurate reading of *SpO₂* values is among the most important factors for successful management of LTOT care. With the **GO₂** fingertip oximeter, engineered with Nonin Medical’s proven PureSAT® *SpO₂* technology (the same technology that’s in our medical grade oximeters), you can always trust the oxygen saturation and heart rate values reported by your patients. With over 22 years of manufacturing, engineering and research experience, we have a long-standing reputation for excellence. In fact, 96% of pulmonary rehab facilities use Nonin.

Partnering with your Patients

Having a fingertip oximeter that is always accurate helps patients stay within their prescribed oxygen saturation and heart rate levels — enabling better management of their care by giving them the information they need to adjust their oxygen to meet their activity levels. With the power and confidence to manage their oxygen, the **GO₂** fingertip oximeter empowers patients to respond in situations where they find themselves not feeling 100%. Visit www.go2nonin.com to learn more about how the **GO₂** can benefit your patients.

“Every COPD patient should have their own oximeter. I treat it just the same as I would having their own blood pressure cuff.”

—Brenda Crowe, RT

Promoting Exercise Compliance

With the reassurance of the **GO₂** fingertip oximeter by their side, patients will be more apt to comply with an exercise regimen.

Empowering Patients to Take Control of their Lives

A dependable and accurate tool like the **GO₂** fingertip oximeter will empower patients with the knowledge and confidence that they can manage this disease and still live a fulfilling life.

A Message for All

Not All Oximeters are Created Equal!

- Nonin welcomes scrutiny of our independent testing results and they are available for all to see! Check out why Nonin is the most trusted name in fingertip oximetry at www.go2nonin.com.
- Anyone can make low perfusion claims, but ONLY Nonin has claims for fingertip pulse oximeters cleared by the FDA. Check out our accuracy specifications at www.go2nonin.com. Demand the same from the imitators.
- Nonin’s quality standards require every **GO₂** be able to withstand 25 drops from one meter onto a hard surface. Try it with a **GO₂** imitator and be prepared to pick up the pieces!
- Home environments can get messy — **GO₂** exceeds IP33 water ingress testing.
- The **GO₂** is NOT a throw-away product. If something happens, it can be repaired and like new in record time. Nonin is a company that will be here tomorrow.
- At Nonin, you’ll speak with a real person right away. The voice you hear represents the pride in our products. We go above and beyond to have satisfied customers.