



PATIENT CARE



BLANKET WARMER FROM CHECK CORPORATION

PROVIDES SOOTHING WARMTH TO HOSPITAL AND CLINIC PATIENTS

Hospitals and clinics around the country are adding the Electric HOTBAG Blanket Warmer to their equipment list to provide soothing warmth to patients. When plugged into any 110v AC outlet using the AC power adapter (sold separately), the Electric HOTBAG Blanket Warmer will keep medical blankets at a constant 160 degrees indefinitely, always ready for use.



- Made of 1000 Denier Nylon with Dupont Thermoloft insulation
- 18" Wide X 18" Long X 6.5" High
- Thermostat Controlled Heating Elements (Top & Bottom)
- Electric HOTBAG Interior Heats to 160-175°
- 4.5 Amps Power Consumption
- Weighs just 3 Lbs.
- Requires 110v AC power supply for use indoors (sold separately) Power supply is UL listed

The Electric HOTBAG Blanket Warmer is made of rugged 1000 Denier nylon and Dupont's Thermoloft insulation for maximum heat retention. Each bag has thermostat controlled heating elements in the top and bottom that can be easily removed for cleaning or laundering of the HOTBAG.

ICE RESCUE March 25th, 2005

"Our crew began warming a hypothermia victim with a blanket from your HOTBAG after she fell through the ice. By the time they reached the hospital her core temperature had been raised from the low 80s to 95.1 degrees F."

"Your HOTBAG saved her life"

Brett Murphy, Director of Operations
Regional EMS - Flint, MI

(News story video can be seen on www.HOTBAG.com)



Made in USA

www.HOTBAG.com

800-WARMS-UP (800-927-6787)

Phone: 248-680-2323 Fax: 248-680-2326

available from:

Emech Medical - 0800 436 324

sales@emech.co.nz

www.emech.co.nz

Another Quality Heating Product From CHECK CORP