


INOGEN ONE^{G3} System



Inogen One[®]
AC Power Supply*
Item # BA-301



Inogen One[®] G3
DC Power Cable*
Item # BA-306



Inogen One[®] G3
Carry Bag*
Item # CA-300



**Only
4.9 pounds!**



Inogen One[®] G3
Single Lithium Ion Battery*
Up to 4 hours run time
Item #BA-300

The Inogen One[®] G3 System (IS-300):

- 45% Smaller
- 30% Lighter
- 27% More Oxygen Per Pound
Compared to the Inogen One G2

*Included with system

**Accessory sold separately



Inogen One[®] G3
Double Lithium Ion Battery**
Up to 8 hours run time
Item #BA-316

Oxygen. Anytime. Anywhere[®]

The new Inogen One[®] G3 delivers the independence of previous Inogen One portable concentrators, but it's smaller, lighter, and makes more oxygen per pound. It's a portable concentrator light enough to be worn care-free by the active oxygen user. That means there are no tanks to refill and no worry about running out of oxygen.

The Inogen One[®] G3 is your one solution, your single solution, for oxygen at home or away. The Inogen One[®] G3 is Oxygen. Anytime. Anywhere[®].

To learn more, call toll-free 1-877-845-2064.

Weight	4.9 lbs. (includes single battery)
Size	Length: 8.75" Width: 3.0" Height: 7.25"
Oxygen Flow	Intelligent Delivery Technology Four flow settings 1, 2, 3, 4
Power Sensing	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): 3.5 to 4 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): 7 to 8 hours Recharge approximately 4-8 hours with AC or DC power
Noise Level	42 Decibels*
Warranty	3 year**
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use, at home or away

*At setting 2

**See Inogen One[®] Limited Warranty Statement for complete warranty information

Manufactured by:



UK Distributor
Pure O2,
Viking House,
71 Princess Road,
Urmston,
Manchester. M41 5ST
Int: 0044 (0)161 747 2617
UK: 0870 712 02 02
Email: info@healthoxygen.com
Website: www.healthoxygen.com

