

# LID LED EXAMINATION LAMPS



## LID 7 x LED Lamp

- 50,000 hour life - Never change a bulb again!
- 7w with 7 light sources
- Excellent light without Heat or U.V.
- Cooled by heatsink - No moving parts, no faults, no noise
- 1000mm neck
- 27,000 Lux at 50cm with an illuminated field of 20 cm.
- 5500° Kelvin giving a bright white light.
- Excellent colour rendering.
- No UV output and no infrared
- Suitable for all mounting brackets



## LID 3 x LED Lamp

- 50,000 hour life - Never change a bulb again!
- Very low power usage - 4.2w with 3 light sources
- Excellent light without Heat or U.V.
- Cooled by heatsink - No moving parts, no faults, no noise
- 900mm neck
- 12,000 Lux at 30cm
- White light similar to the daylight :5500° K
- Excellent colour rendering. No UV output and no infrared
- Suitable for all mounting brackets



# LID HALOGEN EXAMINATION LAMPS



## LID Halogen Lamp 20W

- Powerful and comfortable to use
- Switch on the lamp base
- 9,000 Lux at 50cm
- 5,000hr bulb life
- Temperature color : 3000° K
- 900mm neck
- 23,000 Lux at 50cm
- Suitable for all mounting brackets



## LID Halogen Lamp 50W


- Dimmer and positioning handle
- Switch on the lamp base
- 15,000 Lux at 1 meter
- 5,000hr bulb life
- Temperature color : 3000° K
- 1000mm neck
- 50,000 Lux at 50cm
- Suitable for all mounting brackets



## Mounting Solutions For Halogen & LED Lamps




Mobile Stand  
Adjustable height




Wall Bracket



Desk clamp



Rail clamp



Universal clamp



Heavy base 1.5kg



Wall Mount with Articulating Arms 2  
x 400mm 3 rotation points