



## Maxi HT (Heavy Terrain)

The AMW Champion. This model has been developed specifically to meet the hard driving conditions in NZ. Its huge motor, tyres and batteries combined with large chassis mean this scooter can handle tough terrain, big, heavy or tall riders and is a real long distance tourer. This is definitely the best valued high performance scooter in New Zealand.



### Features

<b>Motor size:</b>	1200 Watt
<b>Battery power:</b>	2 x 70 Amp (12v)
<b>Tyre size:</b>	14"
<b>Dimensions:</b>	1540 x 670mm
<b>Turning radius:</b>	1.65m
<b>Ground clearance:</b>	160mm
<b>Max speed:</b>	12km/hr
<b>Range:</b>	35-40km (User pending)
<b>Carry capacity:</b>	150-200kg (Terrain pending)
<b>Climbing grade:</b>	15° degrees
<b>Construction:</b>	Penny & Giles Controllers

<b>Brakes:</b>	Electronic regenerative
<b>Volts:</b>	24v DC
<b>Suspension:</b>	Full
<b>Seat:</b>	Adjustable front/back, swivel, height, backrest.
<b>Quick Release:</b>	Dash mounted brake release
<b>Controller:</b>	P&G
<b>Package Includes:</b>	Front/back baskets, rain cover, manual, walking stick holder, walking stick, charger, paint repair, toolkit.

Available in:

