



Desk model - Long-lasting mercury sphygmomanometer as a desk model in a rugged **metal case** with a two-tone finish in light gray / black. The metal faceplate with easy-to-read scale, the practical bulb support and the **chromium-plated metal parts** bear witness to stylish workmanship and maximum precision. Also optionally available in a latex-free version (when ordering enter LF for latex-free in front of the article code).

- Maximum error tolerance of +/- 3 mmHg
- 99.99% pure mercury
- Precision glass column
Internal diameter 4.2 mm +/- 0.2 mm
- Latex bulb with chromium plated release valve
- Precision air release valve, wear-free and with fine regulation
- Microfilter for long life
- Cuff with two-tube latex bladder
- Metal faceplate with easy-to-read scale up to 300 mmHg

- Mercury lock
- Special seal against mercury contamination
- Large storage compartment for cuff + latex bulb
- Cleaning device for glass column

- No. 1030** Cotton bandage cuff
- No. 1031** Cotton hook cuff
- No. 1032** Calibrated nylon velcro cuff