

MD700 Table Top Oximeter

Features:

- Suitable for both spot and continuous monitoring of SP02 and Pulse
- Extended 3 metre cable connection between finger sensor and monitor available
- Reusable fingertip sensor (paediatric sensor available)
- 3 hour battery backup with rechargeable battery
- Adjustable alarm parameter settings
- Large easy to read LED display
- Accuracy +/- 2 BPM or +/- 2%
- Weighs just 450grams
- Dimensions 190mm (L) x 148mm (W) x 67mm (H)



MD600 Hand Held Oximeter

Features:

- Suitable for both spot and continuous monitoring of SP02 and Pulse
- Reusable fingertip sensor (paediatric sensor available)
- Adjustable alarm parameter settings
- 24 hour data storage
- Dimensions 150mm (L) x 94mm (W) x 47mm (H)
- Accuracy +/- 2 BPM or +/- 2%



MD300 Fingertip Oximeter

Features:

- Fits on to the end of a finger for instant measurement of SP02 and Pulse
- AAA battery operation with auto off
- LED display
- Weighs only 50 grams including batteries
- Accuracy +/- 2 BPM or +/- 2%



MD300C Childs' Fingertip Oximeter

Features:

- - small built in finger sensor to accommodate paediatric patients from 3 months old
- - fits on to the end of a finger for instant measurement of SP02 and Pulse
- - AAA battery operation with auto off
- - LED display
- - accuracy +/- 2 BPM or +/- 2%

