

WALMAY CRYOSPRAY

available from:

Emech Medical

0800 436 324

www.emech.co.nz

sales@emech.co.nz



Cryosurgery:

Patient benefits at your fingertips for both the Dermatologist & General Practitioner.



The Cryospray is held at the balance point, with the fingers falling naturally on the ergonomic regulator control. This makes it a comfortable and fatigue free instrument to use.

The Advantage of Using Walmay Cryospray

Primary care physicians can simply and quickly treat such problems as keratoses, condylomata, warts, skin tags and other skin lesions every day using **Walmay Cryospray**.

Not only is the procedure simple, there are many benefits for both the physician and the patient in its use.

These include:-

- Minimum preparation and quick application, therefore faster patient turn around
- No anaesthesia
- Little if any patient after care required
- Cryosurgery is well received by patients

Liquid nitrogen cryosurgery using **Walmay Cryospray** is probably the simplest and quickest procedure you could execute during a consultation. Liquid Nitrogen is relatively inexpensive and widely distributed throughout Australia.

In times of increasing cost pressures for both the physician and patient, the practitioner who uses **Walmay Cryospray** benefits both by:-

- Increasing productivity. In about the same time that it takes to write a referral, the primary care practitioner can treat the patient.
- Saving the patient considerable time and expense. The treatment is normally carried out in the one interview, therefore your patient is spared the time and associated expenses of a referral.



Some advanced and complicated conditions you will refer to a dermatologist, but there are many, times when your practice would make use of **Walmay Cryospray**.

By using **Walmay Cryospray** you will be making available another service to your patients not widely offered by primary care physicians. This will build your practice and further help in patient retention.

The **Walmay Cryospray** is designed and manufactured in Australian, precision made instrument ensuring the physician continuing reliability. If you are already a cryosurgery practitioner, you can confidently update with a Walmay Cryospray.

Features of Walmay Cryospray

- The regulator and safety valve, incorporated in the one unit, are our own design and is precision machined from a solid brass billet.
- A stainless steel vacuum container that feels comfortable in the hand with the operating fingers falling naturally on the regulator control. This makes it a balanced and fatigue free instrument to use.
- The 0.3L vacuum container will retain a full days supply of liquid nitrogen.
- Set of five brass spray tips - four directional cone shaped tips, with varying orifice, for accurate visual aiming during treatment, and one bent spray extension are supplied with the instrument.
- Registered under the ARTG Therapeutic Goods Act 1989.
- Table stand for product stability.
- Carry bag for easy transport from surgery to surgery.
- Designed and manufactured in Australia.

Accessories

- Extractor, available for unit filling.

WALMAY HOLDINGS PTY LTD
WALMAY

9 Artisan Road, Seven Hills NSW 2147
PO Box 144, Kellyville NSW 2155

P. 02 9624 3188
F. 02 9624 4387
E. info@walmay.com.au