


AARON® 2250™

**Bovie®**

emech  
+ + +  
m e d i c a l  
0800 436 324 | sales@emech.co.nz

# ELECTROSURGERY FOR THE HOSPITAL OPERATING ROOM AND SURGERY CENTER



ELECTROSURGICAL GENERATOR

## CUTTING EDGE

200 watts of Cutting power

## COMPETENCE

Enough power to perform any procedure in a surgery center and most Hospital OR procedures

## VERSATILITY

Offers Pure Cut, ten (10) modes of Blend, two (2) modes of Coagulation, Bipolar Coagulation and BovieFCFS™ technology allowing the Doctor to move almost effortlessly from Monopolar to Bipolar during procedure

## AFFORDABLE

The highest in safety and quality were built into a 200 watt generator for maximum reliability with economy

## SECURITY

The Aaron 2250 is designed with state-of-the-art safety features including Bovie NEM™ pad sensing technology which monitors the return electrode for optimum patient protection. BovieDED™ monitors all functions of the generator to ensure constant safety for the patient, doctor and OR staff





## RELIABILITY

The Aaron 2250 is backed by Bovie's two (2) year manufacturer's limited warranty

OTHER OPTIONAL  
ACCESSORIES

Disposable & Reusable Pencils*	Disposable Grounding Pads*	Monopolar & Bipolar Footswitches
		
ESP1 & ESPR	ESRE & ESREC	BV-1253B & BV-1254B

Disposable Suction Coagulator*	Reusable SuperCut™ Tungsten Needle*	Resistick II™ Coated Blade*
		
SCH10 Handswitch	ES60R	ES01T

Disposable and Reusable Electrodes*		Bipolar Forceps*	Bipolar Forcep Cord*
Blade	Ball	Bayonet	Reusable
			
ES01	ES20	A821	A827V

\*More styles available

Industry-standard accessories fit most electrosurgical generators

Additional accessories available — contact your Bovie representative for details.

TECHNICAL  
SPECIFICATIONS

Technical Specifications...

<b>Dimensions:</b>		<b>Line Voltage:</b>	<b>Line Frequency:</b>
Height	15.30 cm (6.00 in.)	100-240 VAC	50 - 60 Hz
Width	31.10 cm (12.30 in.)		
Depth	41.30 cm (16.30 in.)	<b>Output Type</b>	<b>Duty Cycle</b>
Weight	<8.8 kg (< 19 lbs)	RF ISOLATED <input checked="" type="checkbox"/> F	10 S/30 S
<b>Monopolar</b>	<b>Output Power</b>	<b>Output Frequency</b>	<b>Rep Rate</b>
Cut	200 W @ 300 Ω	490 kHz ± 5 kHz	N/A
Blend (MAX)	200 W @ 300 Ω	490 kHz ± 5 kHz	30 kHz ± 5 kHz
Coagulation	120 W @ 500 Ω	490 kHz ± 5 kHz	30 kHz ± 5 kHz
Fulguration	80 W @ 500 Ω	490 kHz ± 5 kHz	30 kHz ± 5 kHz
Bipolar	80 W @ 150 Ω	490 kHz ± 5 kHz	30 kHz ± 5 kHz
<b>Warranty:</b> limited two (2) year on unit only			

Aaron® 2250™ Unit



BV-IDS-CS Stand  
(Optional)

Smoke Evacuation System  
(Optional)



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com  
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie® and Aaron® are registered trademarks of Bovie Medical Corporation.