

e-mega

precision meets efficiency

SPHYGMOMANOMETERS



e-mega aneroid sphygmomanometer

State of the art aneroid palm style sphygmomanometer available as 1- and 2-tube version with disinfectable latex-free one-piece cuffs.

- **Fast inflation of the cuff due to new ergonomic spoon design**
 - **Lifelong precision due to almost non-aging precision movement, with specially tempered copper-beryllium diaphragm for pressure-loading capacity up to 600 mm Hg**
 - **Disinfectable and washable one-piece cuffs in different sizes**
 - **Comfort – air release valve, wearfree, with fine adjustment**
- Large manometer scale (ø 63 mm) and optimized scale and pointer design for better readability of values
 - High-polished manometer surface, easy to disinfect
 - Maximum error tolerance +/- 3 mm Hg
 - Microfilter protects air release valve and measuring system
 - Available as 1-tube and 2-tube model with plug-in connection
 - Completely latex-free version with shock-protecting polycarbonate casing
 - Delivery includes one-piece cuff and vinyl zipper bag
- e-mega with disinfectable bladderless cuff for adults, 1-tube
 - e-mega with disinfectable bladderless cuff for adults, 2-tube

white Art. No.: 1370-150 black Art. No.: 1375-150
 white Art. No.: 1371-151 black Art. No.: 1376-151



Rudolf Riester GmbH | P.O. Box 35 | Bruckstr. 31 | D-72417 Jungingen | Germany
 Tel.: (+49) +7477-9270-0 | Fax: (+49) +7477-9270-70 | info@riester.de | www.riester.de

0800 436 324 | www.emech.co.nz | sales@emech.co.nz