


BOWIE & DICK TEST e HELIX TEST process challenge devices



Complying Test.

BOWIE & DICK TEST

	U.M.	
Overall dimensions W x H x D	mm	120 x 28 x 120
Weight (mass)	g	220±4
Room temperature allowed	°C	from 0 to 35
Relative humidity allowed	%	from 20 to 60
Pieces per box	n	20

Technical data are subject to change without notice

HELIX TEST

	U.M.	
Length	mm	1500
Internal diameter	mm	2
Chemical indicators per box	n	100

Technical data are subject to change without notice

BOWIE&DICK TEST – HELIX TEST can check the complete sterilization phase. The security system can be customized: **chemical process indicators** are aimed at verifying the sterilization phase. **Quality** and **reliability** are **guaranteed** by Cominox brand.