

GENERAL SKIN EXAMINATION WITH POLARIZED LIGHT

Lumio is the ideal device for general skin exams where an in-depth, glare-free image of the skin is essential.

With its 75 mm lens offering a magnification of 2x, the power of 40 bright-white LEDs (light emitting diodes), and cross-polarization, Lumio is a great tool for a variety of dermatological exams: Varicose

veins, pigmented skin lesions, hair follicles and much more.

Lumio is powered by four AA batteries that allow for a minimum of 250 exams at 30 seconds each.

Model No. LUM is supplied with one Lumio, four AA batteries and one protective neoprene pouch.



Lumio®
dermLite®



REF: LUM • PERMANENT CROSS-POLARIZATION • 40 LEDS • AA BATTERIES • 75 MM 2X LENS •
STYLISH, ERGONOMIC DESIGN • PROTECTIVE NEOPRENE POUCH