

# DERMLITE **DL1**

## Incredibly versatile digital dermoscopy.

DL1 packs the innovations that make the top-of-the-line DermLite devices so convenient & powerful into a digital dermoscopy solution designed for your smartphone or tablet.

Our designers squeezed a superb, full-size 15 mm lens that's optimized for iPhone, high-powered LED lighting and a high-performance rechargeable battery

into an all-aluminum design smaller than your ring finger. Snap it to your iPhone, Samsung® Galaxy, iPad or iPad mini, and capture images of excellent clarity and sharpness in both polarized and non-polarized modes.

While DL1 ships with an iPhone 6/6S adapter, inexpensive adapters for above devices are sold separately,



iPad mini



and since we are constantly developing other adapters for future smartphones and tablets, upgrades are easy.

Conveniently, its battery can be recharged using any USB port or charger and, due to its unique lighting angle, it works just as well for ENT exams as it does for dermoscopy. That's versatility with vision.



iPad



DL1 & iPhone 6

- + Works without a mobile device
- + Includes adapter for iPhone 6/6S
- + Works on a variety of mobile devices
- + Cross-polarized and non-polarized modes
- + Contact & non-contact dermoscopy
- + Glass faceplate with 10 mm markings
- + Recharges from any USB jack or charger
- + Used without phone, it clips to your shirt

