




Remote Display can be mounted on a wall or placed on a table

Advantages

- EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Remote Display**
The remote display provides versatility with the ability to be mounted on a wall or placed on a table
- 600 lb Capacity Floor Scale**
The 752KL is the market's highest capacity floor scale, providing more versatile use
- Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)
- Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience





PORTROD
Optional Wall Mounted
Height Rod Available



Remote Display can be mounted on a wall or placed on a table

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
752KL	25.625" x 16.75" x 4.625" / 1.15	18 lb	892076002138

Specifications



- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 14 1/4" (w) x 14 1/4" (d) x 2 5/8" (h) / 362 mm (w) x 362 mm (d) x 67 mm (h)
- Height Rod: Optional PORTROD
- Display: 1" LCD
- Connectivity: USB
- Power Source: 110-240V adapter included, 6 C-cell batteries not included
- Optional Carrying Case: 64771
- Functions: LB / KG Conversion, LB / KG Lock Out , Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- 2 Year Limited Warranty