


Health o meter®
P r o f e s s i o n a l

10 Year Limited
Warranty



Advantages

- ✓ **Powder Coated Fabrication**
Powder coating inside and out protects it from rust, extending the life of the scale
- ✓ **Height Rod Included**
The easy-to-use height rod provides height measurements from 23" to 84" / 60 cm - 213 cm
- ✓ **Exclusive 10 Year Limited Warranty**
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship



0800 436 324 | 1800 043 324
www.emechmedical.com
sales@emechmedical.com

Health o meter[®]
P r o f e s s i o n a l



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
402KL	58.875" x 20.25" x 6.375" / 4.40	33.25 lb	892076002015

Optional Wheels
Item # 55000



Optional Counterweights
Item # 55090



Specifications

- Capacity: 390 lb / 180 kg
- Graduation: 1/4 lb / 100 g
- Platform Size (w x d x h): 10 1/2" x 14" x 3 1/4" / 267 mm x 356 mm x 83 mm, platform mat: non-slip black
- Fixed Beam, pounds and kilograms can be read on both sides
- Professional Dual Reading Height Rod: 23 5/8" - 84" / 60 cm - 213 cm
- Optional wheels available (item # 55000)
- Optional counter weights available, to increase capacity to 490lb / 210kg (item # 55090)
- 10 Year Limited Warranty

emech