



Advantages

- 440 lb / 200 kg Capacity Floor Scale**
This high capacity floor scale is versatile for most patients
- 10-Second Weight Hold**
Weight measurement stays displayed for 10 seconds after stepping off the scale so patients can stabilize themselves
- Large 2.5" LCD Display**
The large display offers easy-to-read weight measurement
- Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more stable weighing experience
- Wide Platform**
The spacious platform allows patients to find a comfortable and secure stance to help maintain their balance
- Non-Slip Surface**
The platform's textured mat provides traction and a stable foothold

Health o meter®
P r o f e s s i o n a l



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
880KL	4	18.5" x 14.4" x 6.9" / 1.06	42 lb	854896003665



MADE IN HUNGARY



Specifications

- Capacity: 440 lb / 200 kg / 31 st
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 17" (w) X 13" (d) X 1.3" (h) / 432 mm (w) x 330 mm (d) x 33 mm (h)
- Display: 2-1/2"
- Power Source: (2) 3-volt lithium batteries (included)
- Functions: LB / KG / ST Switch, Auto Off
- 10 Second Weight Hold
- 1 Year Limited Warranty