


**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l



### Advantages

- ✓ **440 lb / 200 kg Capacity Floor Scale**  
The high capacity of this floor scale provides more versatile use
- ✓ **1-1/4" LCD Display**  
The large display shows the user an easy-to-read weight measurement
- ✓ **EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter<sup>®</sup> Professional Connectivity Solutions
- ✓ **Low Profile Platform**  
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience



0800 436 324 | 1800 043 324  
www.emechmedical.com  
sales@emechmedical.com

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l  
www.homscales.com

**Healthometer®**  
Professional



844KL Display



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
844KL	2	15.125" x 14.75" x 6" / 0.77	8.5 lb	10892076002401



**Specifications**

- Capacity: 440 lb / 200 kg
- Resolution: 0.1 lb / 0.05 kg
- Platform Size: 12 5/8" (w) x 12 5/8" (d) x 1 7/8" (h) / 321 mm (w) x 321 mm (d) x 48 mm (h)
- Display: 1 1/4" LCD
- Power Source: (2) 3-Volt Lithium batteries included
- Connectivity: USB
- Optional Carrying Case item # 64771
- Functions: LB / KG Switch, Auto Zero,
- 5 Second Weight Hold, Auto Off
- 1 Year Limited Warranty



Carrying Case for 844KL  
Item# 64771

