

Remote Display can be mounted on  
a wall or placed on a table



## Advantages

✓ **Durable Construction**

The scale has an easy-to-clean stainless steel solid construction and durable, non-skid mat and can withstand frequent use

✓ **600 lb Capacity Floor Scale**

The high capacity of this floor scale provides more versatile use

✓ **AC Adapter Included**

Scale is powered by a 120V AC adaptor (included) or 6 C-Cell batteries (not included) which allows for the convenience of dual power sources

✓ **Motion Detection**

Scale has the ability to weigh accurately despite movement of the subject on the scale

emech

m e d i c a l

0800 436 324 | 1800 043 324  
sales@emechmedical.com  
www.emechmedical.com

Remote Display can be mounted  
on a wall or placed on a table



2842KL Display



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2842KL	1	47.5" x 25.625" x 5.25" / 3.70	61.75 lb	892076002480



### Specifications

- Capacity: 600 lb / 270 kg
- Resolution: 0.2 lb / 0.1 kg
- 7/8" LCD display
- Functions: LB / KG Conversion, Zero Out / Tare, Motion Detection
- Platform Size: 22¼" (w) x 42" (d) x 2½" (h) / 565 mm (w) x 1067 mm (d) x 64 mm (h)
- 120V adapter included, 6 C-cell batteries (not included)
- 3 Year Limited Warranty