

## Class 4 Dual Multi Parameter Indicator

The Valisafe Class 4 Dual Multi Parameter indicators offer an economical yet effective way of monitoring your steam sterilization processes.

(i) With one accurate calibration world wide it is a secure and safe monitoring system yet simple and easy to use.

(ii) Unlike most other class 4 indicators available today the Valisafe indicator is non toxic - another sign of it's advanced technology.

- Conforms to EN ISO 11140–1 Class 4
- Class 4 Multi Parameter indicator - Steam
- Monitors all 121°C to 134°C steam sterilization cycles
- Clear accurate colour change
- Easy to interpret - easy to record
- For use in every tray or item
- Non-hazardous / Lead free according to EC Regulation 07/2006
- Strips dividable to make 500
- Laminated

<b>Name:</b>	Dual Multi Parameter Indicator
<b>Classification:</b>	EN-ISO 11140-1 (Class 4)
<b>Operating temp:</b>	121°C -134°C
<b>Storage conditions:</b>	10-30°C
<b>RH:</b>	30-70°C
<b>Correct symbols:</b>	In accordance with EN ISO 11140-1
<b>Indicator:</b>	Class 4
<b>Disp/single use:</b>	Yes
<b>Coated:</b>	Yes
<b>Non-toxic/Lead free:</b>	Yes, according to EC Regulation 1907/2006
<b>Lot no.</b>	Yes
<b>Prod date/Exp date:</b>	Yes
<b>Operator:</b>	Yes
<b>Pass/fail result:</b>	Yes
<b>Colour change:</b>	Blue to dark grey or black



Before  
Original indicator



After  
After total process



Order details:

V3502038	Dual Multi Parameter Indicator	250 Indicators / can be split (500)
----------	--------------------------------	-------------------------------------