

iMAC 120

High-end 12-Ch ECG machine



Zoncare
FDA FC CE ISO9001

WUHAN ZONCARE BIO-MEDICAL ELECTRONICS CO.,LTD.

Address: Zoncare Building, No. 380, High-tech 2nd Road, Eastlake High-tech
Development Zone, Wuhan, Hubei 430206 P. R. China.

Tel: +86-27-8777 0581 Fax: +86-27-8777 0203

Email: info@zoncare.com <http://www.zoncare.com>

© 2017 Wuhan Zoncare Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.



<http://www.zoncare.com>






Excellent Performance

- Digital Sampling Rate: 32000Hz
- A/D Conversion: 24 bits
- 1s baseline stability
- Frequency Response: 0.01-350Hz
- Support sending ECG report by Email
- Support double-pole pacemaker detection
- Glasgow ECG analysis algorithm
- ST segment elevation myocardial infarction diagnosis ability is in the industry leading position
- Original iFilter™ adaptive filtering technology has greatly improved the ECG signal quality



Recording weak signal accurately




Single-Chamber Pacing



Dual-Chamber Pacing



IMAC ECG Machine




Other Brand ECG Machine

Innovative Ergonomic Design

- 10.4 inch super sensitive and capacitive touch screen design
- Virtual and physical button input
- Innovative handle design , support screen angle
- Standard configured the multi-axis TP® innovative thermal printer can accurately trace each ECG waveform



Drawer type thermal paper recording device



Three-pronged ECG Leads, Effectively reduce winding



Trolley

