

Bovie® Bantam | PRO™

Electrosurgical Generator + High Frequency Desiccator

Technical Specifications

MONOPOLAR	WATTS	VPEAK MAX VOLTAGE	CREST FACTOR (TYPICAL)	FREQUENCY (KHZ)
Cut I	50	1600V	2.2 ± 20%	343kHz ± 10%
Blend	50	2100V	3.5 ± 20%	368 kHz ± 10%
Coagulation	50	2900V	5.2 ± 20%	340 kHz ± 20%
Spray Coagulation	40	6300V	9.5 ± 20%	410 kHz ± 10%

MONOPOLAR	WATTS	VPEAK MAX VOLTAGE	CREST FACTOR (TYPICAL)	FREQUENCY (KHZ)
Bipolar	40	950V	5.5 ± 20%	368kHz ± 10%
Micro Bipolar	40	300V	2.6 ± 20%	338 kHz ± 10%

Line Frequency	50-60 Hz
Output Frequency	304-451 kHz
Line Voltage	100-240 VAC
Cooling System	Passive Cooling
Technology	Isolated RF Output
Monopolar HP Ports	1 Port(s)

Warranty	4 years
Width	228 mm (8.98 in.)
Range Trip Point	(10C° to 40°C)
Depth	105 mm (4.13 in.)
Weight	< 2.26 kg (< 5 lbs.)
Height	188 mm (7.40 in.)
Certifications	CE 0413